




SPIGGLE & THEIS
Medizintechnik

Sinus Instruments



Quality

Made in
Germany

Shaping the future of ENT care
smarter, safer and beyond limits



About Us

WELCOME TO OUR INSTRUMENT CATALOG

Before you immerse yourself in the world of our “MADE IN GERMANY” instruments, allow us to offer a brief glimpse into our success story.

What began as a modest trading company three decades ago has since developed into a globally recognized manufacturer of top-tier ENT instruments and consumables.

Our history is defined by a profound dedication to ENT surgery, healing practices, and continuous growth. Our passion resonates perfectly with the physicians, specialists, and medical professionals who are integral to our network and innovative spirit.

Another vital aspect of our success is our employees, some of whom have been with our company for decades. Their expertise and unwavering dedication have empowered us to collaborate with you in developing and consistently enhancing new products.

In 2020, SPIGGLE & THEIS Medizintechnik GmbH was primarily acquired by Thyssen'sche Handelsgesellschaft, ensuring the company's continued expansion while prioritizing our innovative capacity and commitment to quality. With this robust partnership at our side, new horizons have emerged, accelerating the realization of our visions and enhancing our ability to meet customer demands more efficiently.

This catalog showcases premium instruments tailored to the needs of our clients and partners in the field of ENT medicine. We hope you discover the perfect instrument that meets your requirements and allows our commitment to quality to speak for itself.

Thank you very much for your loyalty and trust.

Your SPIGGLE & THEIS Team



For hearing, breathing and voice

ENT care for life



Glossary

TL	Total Length	OD	Outer Diameter
WL	Working Length	ID	Inner Diameter
TH	Total Height	Ø	Diameter
TD	Total Depth	WE	Working End
Dim.	Dimension		

Attention

Not every instrument is available in all markets. Availability depends on local regulations. Please contact your SPIGGLE & THEIS representative or distributor, if you have any questions about the availability of our products in your region.

Specula.....	7
Retractors	10
Probes.....	12
Suction & Irrigation.....	14
Knives.....	17
Elevators & Raspatories.....	19
Chisels	22
Hammers	23
Curettes	24
Dissector	28
Scissors	29
Nasal Forceps	33
Grasping Forceps	44
Forceps.....	51
Punches	52
Micro Instruments	60
Outpatient Nose Instruments.....	61
Miscellaneous.....	62




Hartmann nasal speculum

dismountable

Article No.

20-100-01	for adults
20-100-02	for juveniles & children

**Killian-Struycken septum speculum**

Total Length: 14.0 cm • dismountable

Article No.**Working End**

20-150-40	45.0 mm length
20-150-55	55.0 mm length
20-150-65	65.0 mm length
20-150-75	75.0 mm length
20-150-90	90.0 mm length



Specula

Sinus speculum

for children • Total Length: 13.0 cm
Working End: left blade 25 mm • dismountable
Article No.: 30-702-25



Sinus speculum


for children • Total Length: 14.5 cm
Working End: right blade 25 mm • dismountable
Article No.: 30-703-25



Sinus speculum

for adults • Total Length: 13.5 cm • dismountable

Article No.	Working End
30-705-55	left blade 55 mm
30-705-65	left blade 65 mm
30-705-75	left blade 75 mm
30-706-55	right blade 55 mm
30-706-65	right blade 65 mm
30-706-75	right blade 75 mm



Selfholding Spekulum

with adjusting screw • Total Length: 15.0 cm
Working End: left and right slotted • dismountable
Article No.: 30-711-00



Flexible holding arm for use

with selfretaining nasal specula • incl. holding rod

Article No.: 30-700-00



Rotation clamp

for fixation at operating table,
compatible with REF 30-700-00 & 40-337-00

Article No.: 30-700-01



Retractors

Cottle alar retractor

Article No.	TL	Working End
20-560-14	14.0 cm	short blade
20-560-15	16.5 cm	
20-561-15	16.0 cm	long blade



Kocher-Langenbeck retractor

ring handle


Article No.	TL	Dim.
80-308-07	21.0 cm	35.0 × 7.0 mm
80-308-10	22.0 cm	36.0 × 10.0 mm
80-308-15	22.0 cm	35.0 × 15.0 mm
80-308-20	23.0 cm	61.0 × 20.0 mm



Langenbeck-Green retractor

Total Length: 17.0 cm • Dimension: 25.0 × 6.0 mm • flat handle

Article No.: 80-334-16



Middeldorpf retractor

ring handle


Article No.	TL	Dim.
80-306-17	22.0 cm	15.0 × 16.5 mm
80-306-22	23.0 cm	20.0 × 22.0 mm
80-306-30	24.0 cm	25.0 × 30.0 mm



Micro hook

Total Length: 21.5 cm • Working End: straight, blunt forward / backward


Article No.: 90-214-00



Fisch retractor

with joints • Total Length: 16.5 cm • Working End: semi-sharp 3:3 prongs

Article No.: 10-468-17




Probes

Enucleator

Total Length: 27.0 cm • Shaft: bayonet
Working End: left curved • knurled handle


Article No.: 90-204-01



Nasal probe

Total Length: 17.5 cm • Shape: s-shaped • bendable

Article No.: 20-313-00



Ritter frontal sinus bougie

Total Length: 15.5 cm


Article No.	Ø
30-121-25	2.5 mm
30-121-30	3.0 mm
30-121-40	4.0 mm
30-121-50	5.0 mm



Frontal sinus seeker

Working End: double ended


Article No.	TL	Working End
30-207-03	22.0 cm	90° strongly curved
30-207-01	20.5 cm	slightly curved



Probe hook

Total Length: 27.0 cm • Shaft: bayonet • Working End: blunt

Article No.	Working End
90-208-01	left curved
90-208-02	right curved



Maxillary sinus seeker

Total Length: 19.0 cm • Working End: double ended

Article No.: 30-207-02




Ball-end probe

Total Length: 24.0 cm • Shaft: bayonet

Working End: 9.5 mm length • probe angled

Article No.: 90-213-00




Suction & Irrigation

Suction- / Irrigation tube

Working Length: 17.0 cm • OD: 2.5 mm / 2.0 mm

Article No.: 30-749-15



Suction- and irrigation sheath


Article No.		TL	Working End
30-164-00	for endoscope 4 mm	16.0 cm	0°
30-167-00	for endoscope 2.7 mm	17.5 cm	0°
30-167-30	for endoscope 2.7 mm	17.5 cm	30°



Eicken antrum cannula

Total Length: 14.0 cm • Shaft: long curve


Article No.	OD	Description
30-100-25	2.5 mm	bendable
30-100-30	3.0 mm	bendable
30-100-40	4.0 mm	bendable
30-101-25	2.5 mm	
30-101-30	3.0 mm	
30-101-40	4.0 mm	



Eicken antrum cannula

Total Length: 14.0 cm • Shaft: short curve

Article No.	OD
30-102-25	2.5 mm
30-102-30	3.0 mm
30-102-40	4.0 mm



Sinus suction tube

Total Length: 16.5 cm • bendable


Article No.	OD
30-106-25	2.5 mm
30-106-30	3.0 mm
30-106-40	4.0 mm



Ferguson suction tube

with suction control


Article No.	WL	OD
30-081-06	9.0 cm	6.0 FR
30-081-08	9.0 cm	8.0 FR
30-081-09	10.0 cm	9.0 FR
30-081-10	10.0 cm	10.0 FR
30-081-12	10.0 cm	12.0 FR
30-081-15	10.0 cm	15.0 FR



Frazier suction tube

with suction control • Working Length: 9.0 cm


Article No.	OD
30-080-05	5.0 FR ○
30-080-06	6.0 FR ○
30-080-07	7.0 FR ○
30-080-08	8.0 FR ○
30-080-09	9.0 FR ○
30-080-10	10.0 FR ○
30-080-12	12.0 FR ○
30-080-15	15.0 FR ○



Frazier suction tube

Working Length: 10.0 cm • Working End: angled, ball ended


Article No.	OD
30-082-08	8.0 FR
30-082-10	10.0 FR
30-082-12	12.0 FR
30-082-15	15.0 FR



Frazier suction tube

with suction control • Working Length: 15.0 cm

Article No.	OD
30-747-06	6.0 FR
30-747-07	7.0 FR



Suction & Irrigation

Frazier suction tube

Working Length: 8.0 cm • Shaft: conical
tear formed suction control


Article No.	OD
30-105-06	6.0 FR ☐
30-105-08	8.0 FR ☒
30-105-10	10.0 FR ☐



Frazier suction tube

with suction control • Working Length: 15.0 cm


Article No.	OD
30-747-08	8.0 FR
30-747-09	9.0 FR
30-747-10	10.0 FR
30-747-12	12.0 FR
30-747-15	15.0 FR



Frazier suction tube

with suction control • Working Length: 10.0 cm • OD: 10.0 FR


Article No.	Working End
30-084-01	right curved
30-084-02	left curved



Frazier suction tube

with suction control • Working Length: 10.0 cm • OD: 15.0 FR

Article No.: 30-103-15



Cottle septum knife

Total Length: 14.0 cm


Article No.: 20-500-13



Draf flap knife

Total Length: 19.0 cm


Article No.: 30-744-19



Draf round knife

Total Length: 19.0 cm • Working End: 1.5 mm width

Article No.: 30-743-19



Masing pointed knife

Total Length: 14.0 cm

Article No.	Working End
20-501-14	
20-502-14	smooth




Separator

Total Length: 27.0 cm • Shaft: bayonet

Working End: curved upwards

Article No.	Working End
90-210-09	1.5 mm width
90-210-10	3.5 mm width



Knives

Draf sickle knife

Total Length: 19.0 cm • Working End: curved


Article No.: 30-742-19



Sickle knife

Total Length: 19.5 cm


Article No.	Working End
30-200-01	pointed
30-200-02	rounded
30-200-03	curved, left, pointed
30-200-04	right curved, pointed



Sickle knife

Total Length: 27.0 cm • Shaft: bayonet • Working End: pointed

Article No.	Working End
90-202-01	right curved
90-202-02	left curved



Freer elevator double ended

Working End: sharp-blunt


Article No.	TL
20-306-20	20.0 cm
20-305-18	18.5 cm



Gourney suction elevator

Total Length: 19.0 cm • Working End: 3.0 mm width

Article No.: 20-320-19



Howarth elevator double ended

Total Length: 21.0 cm • Working End: sharp-blunt


Article No.: 20-315-21



Masing elevator/raspatory

Total Length: 21.5 cm

Article No.: 20-300-22



Weber elevator double ended

Total Length: 18.5 cm • Working End: sharp-blunt

Article No.: 20-307-18




Elevators & Raspatories

Draf micro raspatory

Total Length: 25.0 cm

Article No.	Working End
30-740-10	slightly curved, 1.0 mm width
30-740-30	slightly curved, 3.0 mm width
30-740-31	strongly curved, 3.0 mm width
30-740-40	angled, 4.0 mm width



Raspatory

Total Length: 27.0 cm • Shaft: bayonet • Working End: sharp

Article No.	Working End
90-206-01	left curved
90-206-02	right curved
90-206-03	curved downwards
90-206-04	curved upwards



Frontal sinus suction raspatory

with suction control • Working End: curved

Handle length 10.0 cm • bendable

Article No.: 30-108-00



Tönnis raspatory

Total Length: 24.0 cm • Working End: 3.0 mm / 5.0 mm width, double ended


Article No.: 30-759-21



Guillen suction raspatory

Total Length: 15.0 cm • 4.0 mm

Article No.	Working End
20-323-01	straight
20-323-02	left curved
20-323-03	right curved



Guillen suction raspatory

Total Length: 21.5 cm • 4.0 mm

Article No.	Working End
20-322-01	straight
20-322-02	left curved
20-322-03	right curved




Chisels

Cottle chisel

Total Length: 18.5 cm • Shape: flat

Article No.	Working End
20-405-03	3.0 mm width
20-405-04	4.0 mm width
20-405-06	6.0 mm width
20-405-07	7.0 mm width
20-405-09	9.0 mm width
20-405-10	10.0 mm width
20-405-12	12.0 mm width



West chisel

Total Length: 18.0 cm • Working End: 4.5 mm width

Article No.: 30-817-40



West gouge

Total Length: 18.0 cm • Working End: 4.0 mm width


Article No.: 30-816-40



Cottle mallet

Total Length: 18.0 cm

Article No.: 20-400-18



Hajek metal mallet

Total Length: 20.0 cm • 200 g

Article No.: 20-402-21




Curettes

Antrum curette

Total Length: 19.0 cm

Article No.	Description
30-203-01	forward cutting
30-203-02	backward cutting



Antrum curette

Total Length: 19.5 cm • Dimension: 7.0 × 10.5 mm
Shaft: curved • Working End: oval

Article No.: 30-205-01



Antrum curette

Total Length: 19.5 cm • Working End: oblong

Article No.	Dim.	Description
30-205-02	2.5 × 6.0 mm	small
30-205-03	4.0 × 9.0 mm	large



Halle antrum curette

Total Length: 21.0 cm • Working End: double ended

Article No.: 30-805-00



Halle antrum curette

Total Length: 21.0 cm • bendable


Article No.	Size
30-807-01	1
30-807-02	2
30-807-03	3



Spoon curette


Total Length: 27.0 cm • Shaft: bayonet • Working End: straight, upwards, semi-sharp, 2.0 mm width

Article No.: 90-201-01

**Spoon curette**


with 5 mm bar • Total Length: 27.0 cm • Shaft: bayonet
Working End: curved upwards, 90°, semi-sharp

Article No.	Working End
90-201-02	2.0 mm width
90-201-03	3.0 mm width

**Spoon curette**


with 5 mm bar • Total Length: 27.0 cm • Shaft: bayonet
Working End: semi-sharp

Article No.	Working End
90-201-04	curved downwards, 45°, 2.0 mm width
90-201-05	curved upwards, 45°, 2.0 mm width
90-201-06	curved downwards, 45°, 4.0 mm width
90-201-07	curved upwards, 45°, 4.0 mm width

**Spoon curette**

Total Length: 27.0 cm • Shaft: bayonet
Working End: slightly curved, 45°, upwards, semi-sharp

Article No.	Working End
90-201-08	2.0 mm width
90-201-09	4.0 mm width




Curettes

Spoon curette

with 5 mm bar • Total Length: 28.5 cm • Shaft: bayonet

Working End: semi-sharp

Article No.	Working End
90-201-14	curved downwards, 45°, 2.0 mm width
90-201-15	curved upwards, 45°, 2.0 mm width
90-201-16	curved downwards, 45°, 4.0 mm width
90-201-17	curved upwards, 45°, 4.0 mm width




Spoon curette

Total Length: 30.0 cm • Shaft: bayonet

Working End: slightly curved, 45°, semi-sharp


Article No.	Working End
90-201-18	upwards, 2.0 mm width
90-201-19	upwards, 4.0 mm width
90-201-28	curved downwards, 2.0 mm width
90-201-29	curved downwards, 4.0 mm width



Spoon curette

Total Length: 27.0 cm • Shaft: bayonet • Working End: 45°, curved downwards, semi-sharp, 2.5 mm width

Article No.: 90-211-01




Spoon curette

Total Length: 27.0 cm • Shaft: bayonet

Working End: semi-sharp, 4.5 mm width

Article No.	Working End
90-212-01	slightly curved upwards
90-212-02	strongly curved upwards
90-212-03	45°, curved upwards



Ring curette

Total Length: 27.0 cm • Shaft: bayonet • Working End: sharp

Article No.	OD	Working End
90-210-04	2.0 mm	45° curved upwards
90-210-05	2.0 mm	90° curved upwards
90-210-06	4.0 mm	45° curved upwards
90-210-07	4.0 mm	90° curved upwards
90-210-08	5.0 mm	90° curved upwards



Sinus frontal curette

Article No.	TL	Working End
30-205-04	19.0 cm	55° curved, oblong
30-205-06	18.0 cm	90° curved, oblong
30-205-70	18.5 cm	70° curved, oblong



Uffenorde antrum curette

Total Length: 21.0 cm

Article No.	Working End
30-809-01	left curved
30-809-02	right curved




Dissector

Dissector

Total Length: 27.0 cm • Shaft: bayonet

Working End: curved upwards, blunt, 3.0 mm width

Article No.: 90-210-03




Mucosa dissector

Total Length: 24.0 cm • Dimension: 30.0 × 6.0 mm

Shaft: bayonet

Article No.: 90-215-00




Heymann nasal scissors

Premium Line

Premium Line • Working End: straight
serrated

Article No.	TL	Description
20-650-16	16.0 cm	delicate
20-651-18	18.0 cm	robust



20-650-16

Heymann nasal scissors

Working End: straight • serrated

Article No.	TL	Description
20-680-16	16.0 cm	delicate
20-683-18	18.0 cm	robust



Heymann nasal scissors

Total Length: 16.0 cm • serrated

Article No.	Working End
20-681-16	left curved
20-682-16	right curved



20-682-16



Scissors

Walter scissors

Total Length: 10.0 cm • Working End: angled

Article No.: 20-630-10




Conchotomy scissors

Total Length: 11.5 cm • Working End: straight




Article No.	Working End	Description
20-693-00		standard
20-693-01		extra delicate
20-693-02	pointed	extra delicate, lower blade longer



Weber nasal scissors

Total Length: 10.5 cm • Working End: straight
irrigation channel • serrated • golden handle • extra delicate


Article No.: 20-691-01



Draf ethmoid scissors

Dimension: 5.0 × 1.3 × 2.0 mm

Article No.	WL	Working End
30-832-13	13.0 cm	straight
30-832-16	16.0 cm	straight
30-832-19	19.0 cm	straight
30-833-16	16.0 cm	left curved
30-834-16	16.0 cm	right curved
30-835-19	19.0 cm	50° upcurved





Dr. med. Hohenhorst,
Chief physician of the ENT department
of the Alfried Krupp Clinic, Essen,
Germany


“These scissors are ideally suited for opening ethmoidal cells while keeping the mucous membrane largely intact.

The super-fine and elongated lower blade ensures that the mucous membrane of the ethmoid bone is “punctured” and not “pushed forward.”

Nasal concha scissors

Total Length: 13.0 cm • irrigation channel • serrated

Article No.	Working End
20-690-01	straight
20-690-02	left curved
20-690-03	right curved





Advantages Comfort Line

- Quick and easy opening of the sliding shaft
- Validated cleaning
- Safe and comfortable storage
- Made to last

Blakesley nasal forceps

comfort • Working End: straight • golden button

Article No.	WL	Working End
20-915-00	12.0 cm	3.0 mm width
20-915-01	12.0 cm	3.5 mm width
20-915-02	12.0 cm	4.0 mm width
20-915-03	12.0 cm	5.0 mm width
20-915-04	12.0 cm	6.0 mm width
20-915-10	15.0 cm	3.0 mm width
20-915-12	15.0 cm	4.0 mm width



20-915-02



20-915-03



20-915-04



20-915-10

Nasal Forceps

Blakesley nasal forceps

comfort • Working Length: 15.0 cm • Working End: 45° angled, 3.5 mm width • golden button


Article No.: 20-925-11



Blakesley nasal forceps

comfort • Working Length: 12.0 cm • golden button


Article No.	Working End
20-925-00	45° angled, 3.0 mm width
20-925-01	45° angled, 3.5 mm width
20-925-02	45° angled, 4.0 mm width
20-925-03	45° angled, 5.0 mm width
20-925-04	45° angled, 6.0 mm width
20-935-01	90° angled, 3.5 mm width
20-935-02	90° angled, 4.0 mm width



Blakesley nasal forceps

comfort • Working Length: 12.0 cm • Working End: straight throughcutting • golden button


Article No.	Dim.
30-943-01	1.9 × 5.0 mm
30-943-02	2.5 × 6.0 mm
30-943-03	3.0 × 7.0 mm
30-943-04	3.8 × 10.0 mm



Blakesley nasal forceps

comfort • Working Length: 12.0 cm • Working End: 45° angled throughcutting • golden button

Article No.	Dim.
30-943-11	1.9 × 5.0 mm
30-943-12	2.5 × 6.0 mm
30-943-13	3.0 × 7.0 mm



Takahashi nasal forceps

comfort • Working Length: 12.0 cm • golden button

Article No.	Working End
20-945-02	45° curved upwards, 2.5 mm width
20-946-01	straight, 4.0 mm width



Takahashi nasal forceps

comfort • Working Length: 12.0 cm • golden button

Article No.	Working End
20-945-01	straight, 2.5 mm width
20-946-02	45° curved upwards, 4.0 mm width



Nasal Forceps


Advantages

- slim shank
- slender jaw
- smooth finish

Blakesley nasal forceps

Working Length: 11.5 cm • Working End: straight • golden handle


Article No.	Working End
20-215-00	3.0 mm width 
20-215-01	3.5 mm width 
20-215-02	4.0 mm width 
20-215-03	5.0 mm width 
20-215-04	6.0 mm width 



Blakesley nasal forceps

Working Length: 15.0 cm • Working End: straight • golden handle


Article No.	Working End
20-215-10	3.0 mm width 
20-215-12	4.0 mm width 
20-215-14	6.0 mm width 



Blakesley nasal forceps

Working Length: 11.5 cm • Working End: 90° curved upwards
golden handle

Article No.	Working End
20-235-01	3.5 mm width 
20-235-02	4.0 mm width 
20-235-03	5.0 mm width 



Blakesley nasal forceps

Working Length: 11.5 cm • Working End: 45° curved upwards
golden handle

Article No.	Working End
20-225-00	3.0 mm width 
20-225-01	3.5 mm width 
20-225-02	4.0 mm width 
20-225-03	5.0 mm width 
20-225-04	6.0 mm width 



Blakesley nasal forceps

Working Length: 15.0 cm • Working End: 45° curved upwards
golden handle

Article No.	Working End
20-225-11	3.5 mm width 
20-225-12	4.0 mm width 
20-225-14	6.0 mm width 




Nasal Forceps

Blakesley nasal forceps

Working Length: 11.5 cm • Working End: straight throughcutting

Article No.	Dim.
30-340-01	5.0 × 2.0 mm
30-340-02	6.0 × 2.5 mm
30-340-03	7.0 × 3.0 mm



Blakesley nasal forceps

Working Length: 11.0 cm • Working End: straight throughcutting

Article No.	Size	Dim.
30-343-01	1	6.0 × 3.0 mm
30-343-02	2	8.0 × 3.5 mm
30-343-03	3	8.0 × 4.0 mm
30-343-04	4	10.0 × 5.0 mm



Blakesley nasal forceps

Working Length: 11.5 cm • Working End: 45° curved upwards throughcutting


Article No.	Dim.
30-342-01	5.0 × 2.0 mm
30-342-02	6.0 × 2.5 mm
30-342-03	8.0 × 3.0 mm



Blakesley nasal forceps

Working Length: 11.5 cm • Working End: 45° curved upwards throughcutting

Article No.	Size	Dim.
30-343-11	1	8.0 × 3.0 mm
30-343-12	2	8.0 × 3.5 mm
30-343-13	3	9.0 × 4.0 mm
30-343-14	4	10.0 × 5.0 mm



Herberhold cranial through-cutting forceps

Shape: flat

Advantages

- punched material is ejected upwards
- Serration prevents slipping



Article No.	WL	Dim.	Working End
30-332-01	14.0 cm	6.0 × 2.0 mm	straight
30-332-02	14.0 cm	11.5 × 2.5 mm	straight
30-333-01		6.0 × 2.0 mm	60° curved upwards
30-333-02		11.5 × 2.5 mm	60° curved upwards



Nasal Forceps

Through-cutting nasal forceps




Article No.	WL	Dim.	Working End
30-320-01	14.0 cm	5.0 × 1.5 mm	straight 
30-320-02	11.5 cm	6.0 × 2.5 mm	straight 
30-321-01	14.0 cm	5.0 × 1.5 mm	45° angled 
30-321-02	11.5 cm	6.5 × 2.5 mm	45° angled 
30-322-01	14.0 cm	5.0 × 1.5 mm	curved upwards 
30-322-02	14.0 cm	8.0 × 3.0 mm	curved upwards 



Wigand transecting forceps

Working Length: 14.0 cm • Dimension: 4.0 × 3.0 mm • Shape: flat


Article No.	Working End
30-330-01	straight, rounded
30-330-02	45° angled, rounded



Blakesley nasal forceps

Working Length: 13.0 cm • irrigation channel • throughcutting

Article No.	Working End
30-358-00	straight, 3.0 mm width
30-358-01	straight, 3.5 mm width
30-358-02	straight, 4.0 mm width
30-360-00	45° curved upwards, 3.0 mm width
30-360-01	45° curved upwards, 3.5 mm width
30-360-02	45° curved upwards, 4.0 mm width



Through-cutting suction forceps

Working Length: 10.0 cm

Article No.	Shaft	Angle
30-344-00	gerade	
30-344-45		30° angled



Blakesley nasal forceps

Working Length: 16.0 cm


Article No.	Working End
20-240-11	left curved
20-240-12	right curved



Blakesley nasal forceps

Working Length: 11.5 cm

Article No.	Working End
20-240-01	left curved
20-240-02	right curved



Nasal Forceps

Wigand polyp grasping forceps

Working Length: 11.5 cm • Working End: 60° curved upwards

Article No.: 20-236-00



Struempel-Voss nasal forceps

Working Length: 11.5 cm • Working End: 45° curved upwards

Article No.: 20-238-00



Takahashi nasal forceps


Article No.	WL	Working End
20-245-02	11.5 cm	30° curved upwards
20-245-11	15.0 cm	straight, oval
20-245-12	15.0 cm	30° curved upwards, oval
20-246-02	11.5 cm	30° curved upwards, 4.0 mm width



Takahashi suction forceps

Working Length: 10.0 cm • suction tube


Article No.	Working End
30-346-00	straight
30-346-30	30° angled



Blakesley nasal forceps

Working Length: 13.0 cm • irrigation channel

Article No.	Working End
30-350-00	straight, 3.0 mm width
30-350-01	straight, 3.5 mm width
30-350-02	straight, 4.0 mm width
30-350-03	straight, 5.0 mm width
30-352-00	45° curved upwards, 3.0 mm width
30-352-01	45° curved upwards, 3.5 mm width
30-352-02	45° curved upwards, 4.0 mm width
30-354-01	90° curved upwards, 3.5 mm width
30-354-02	90° curved upwards, 4.0 mm width



Grasping Forceps

Grünwald tubular forceps

Working Length: 13.0 cm • Shaft: straight • irrigation channel throughcutting


Article No.	Working End
30-370-15	1.5 mm width
30-370-20	2.0 mm width
30-370-25	2.5 mm width
30-370-30	3.0 mm width



Grünwald tubular forceps

Working Length: 13.0 cm • Shaft: straight • Working End: 45° irrigation channel • throughcutting

Article No.	Working End
30-371-15	1.5 mm width
30-371-20	2.0 mm width
30-371-25	2.5 mm width
30-371-30	3.0 mm width



Grünwald tubular forceps

Working Length: 13.0 cm • Shaft: 30° curved upwards
Working End: straight • irrigation channel • throughcutting

Article No.	Working End
30-372-15	1.5 mm width
30-372-20	2.0 mm width
30-372-25	2.5 mm width




Grünwald tubular forceps

Working Length: 13.0 cm • Shaft: 30° curved upwards

Working End: 45° • irrigation channel • throughcutting

Article No.	Working End
30-373-15	1.5 mm width
30-373-20	2.0 mm width
30-373-25	2.5 mm width




Grünwald tubular forceps

Working Length: 13.0 cm • Shaft: 70° curved upwards

Working End: 2.0 mm width • irrigation channel • throughcutting

Article No.	Working End
30-374-20	70°, straight opening
30-375-20	side to right opening
30-376-20	side to left opening



Frontal sinus forceps

Working Length: 13.0 cm • Working End: 55°, oval cups, 2.0 mm width, horizontal opening • irrigation channel

Article No.: 30-404-02




Grasping Forceps

Frontal sinus forceps

Working Length: 13.0 cm • irrigation channel


Article No.	Working End
30-404-01	55°, oval cups, 2.0 mm width, vertical opening
30-405-01	70°, oval cups, 4.0 mm width, vertical opening



Frontal sinus forceps

Working Length: 13.0 cm • Working End: 2.0 mm width irrigation channel


Article No.	Working End	Description
30-401-01	55°, vertical opening	throughcutting
30-401-02	55°, horizontal opening	throughcutting
30-402-01	90°, vertical opening	throughcutting
30-402-02	90°, horizontal opening	throughcutting



Binner polyp forceps

Working Length: 13.0 cm • Working End: 3.0 mm width, horizontal opening • irrigation channel


Article No.	Working End
30-415-01	45°
30-415-02	80°



Antrum polyp grasping forceps

Total Length: 10.0 cm • Working End: retrograde opening
irrigation channel

Article No.: 30-310-01



Antrum polyp grasping forceps

Total Length: 10.0 cm • Shaft: with extra long shaft curve
Working End: retrograde opening • irrigation channel


Article No.: 30-310-02



Antrum polyp grasping forceps

Total Length: 10.0 cm • Working End: retrograde opening
irrigation channel • rotatable

Article No.	Shaft
30-311-01	
30-311-02	mit extra langem Schaftbogen



Grasping Forceps

Cup forceps

Working Length: 14.0 cm • Dimension: 7.0 × 3.0 mm

Shaft: curved upwards • robust

Article No.: 30-328-01



Flexible alligator forceps

for flexible endoscope 30-171-40 • Working Length: 12.0 cm

Ø: 1.0 mm • Working End: serrated jaws

Article No.: 30-173-40




Flexible biopsy forceps

for flexible endoscope 30-171-40 • Working Length: 12.0 cm

Ø: 1.0 mm • Working End: oval cups

Article No.: 30-172-40



Biopsy forceps

Working Length: 15.0 cm • OD: 2.5 mm • Working End: straight, 4.0 mm length, 2.0 mm width • irrigation channel double-action


Article No.: 30-125-00



Frontal sinus biopsy forceps

Working Length: 13.0 cm • Working End: 3.0 mm width irrigation channel


Article No.	Working End
30-410-01	70°, vertical opening
30-410-02	70°, horizontal opening
30-411-01	110°, vertical opening
30-411-02	110°, horizontal opening



Draf ethmoid forceps

Dimension: 6.0 × 1.0 × 2.0 mm • Working End: straight, finely serrated


Article No.	WL
30-828-16	16.0 cm
30-828-19	19.0 cm



Grasping Forceps

Draf ethmoid forceps


Article No.	WL	Dim.	Working End
30-820-10	10.0 cm	6.0 × 2.5 mm	straight
30-820-13	13.0 cm	6.0 × 2.5 mm	straight
30-820-16	16.0 cm	6.0 × 2.5 mm	straight
30-820-19	19.0 cm	6.0 × 2.5 mm	straight
30-822-13	13.0 cm	6.0 × 2.5 mm	55° angled
30-822-16	16.0 cm	6.0 × 2.5 mm	55° angled
30-822-19	19.0 cm	6.0 × 2.5 mm	55° angled
30-824-13	13.0 cm	6.0 × 2.5 mm	90° angled
30-824-16	16.0 cm	6.0 × 2.5 mm	90° angled
30-824-19	19.0 cm	6.0 × 2.5 mm	90° angled
30-826-16	16.0 cm	4.0 × 1.8 mm	straight



Frontal sinus forceps

Working Length: 13.0 cm • irrigation channel


Article No.	Working End
30-405-02	70°, oval cups, 4.0 mm width, horizontal opening
30-405-11	55°, oval cups, 4.0 mm width, vertical opening
30-405-12	55°, oval cups, 4.0 mm width, horizontal opening
30-407-01	90°, oval cups, 2.0 mm width, vertical opening
30-407-02	90°, oval cups, 2.0 mm width, horizontal opening
30-408-01	110°, oval cups, 4.0 mm width, vertical opening
30-408-02	110°, oval cups, 4.0 mm width, horizontal opening



Cottle columella forceps

Total Length: 11.0 cm

Article No.: 20-700-10




Bipolar forceps


Total Length: 23.0 cm • Shaft: straight, bayonet

Working End: 6.0 mm length, 0.7 mm width • fine handpiece with bar • NON-STICK

Article No.: 80-126-23



Punches



Advantages

- Little effort required
- Smooth running and precise
- ergonomic design
- Totally detachable

Kerrison sphenoid punch

Working Length: 18.0 cm • Working End: 40° curved downwards
totally detachable • golden tip

Article No.	Working End
30-473-10	1.0 mm width
30-473-20	2.0 mm width
30-473-30	3.0 mm width

Kerrison sphenoid punch

Working Length: 18.0 cm • Working End: 40° curved upwards
totally detachable • golden tip • small cutting end

Article No.	Working End
30-469-11	1.0 mm width
30-469-21	2.0 mm width
30-469-31	3.0 mm width
30-469-41	4.0 mm width
30-469-51	5.0 mm width

Kerrison sphenoid punch

Working End: 90° curved upwards • totally detachable
golden tip

Article No.	WL	Working End
30-468-10	18.0 cm	1.0 mm width
30-468-20	18.0 cm	2.0 mm width
30-468-30	18.0 cm	3.0 mm width
30-468-40	18.0 cm	4.0 mm width
30-468-50	18.0 cm	5.0 mm width
30-478-10	23.0 cm	1.0 mm width
30-478-20	23.0 cm	2.0 mm width
30-478-30	23.0 cm	3.0 mm width
30-478-40	23.0 cm	4.0 mm width
30-478-50	23.0 cm	5.0 mm width
30-478-60	23.0 cm	6.0 mm width
30-484-10	20.0 cm	1.0 mm width
30-484-20	20.0 cm	2.0 mm width
30-484-30	20.0 cm	3.0 mm width
30-484-40	20.0 cm	4.0 mm width
30-484-50	20.0 cm	5.0 mm width
30-484-60	20.0 cm	6.0 mm width

Kerrison sphenoid punch

Working End: 90° curved downwards • totally detachable
golden tip

Article No.	WL	Working End
30-471-10	18.0 cm	1.0 mm width
30-471-20	18.0 cm	2.0 mm width
30-471-30	18.0 cm	3.0 mm width
30-479-10	23.0 cm	1.0 mm width
30-479-20	23.0 cm	2.0 mm width
30-479-30	23.0 cm	3.0 mm width
30-479-40	23.0 cm	4.0 mm width
30-479-50	23.0 cm	5.0 mm width
30-479-60	23.0 cm	6.0 mm width
30-485-10	20.0 cm	1.0 mm width
30-485-20	20.0 cm	2.0 mm width
30-485-30	20.0 cm	3.0 mm width
30-485-40	20.0 cm	4.0 mm width
30-485-50	20.0 cm	5.0 mm width
30-485-60	20.0 cm	6.0 mm width

Kerrison sphenoid punch

Working End: 40° curved upwards • totally detachable
golden tip

Article No.	WL	Working End
30-469-10	18.0 cm	1.0 mm width
30-469-20	18.0 cm	2.0 mm width
30-469-30	18.0 cm	3.0 mm width
30-469-40	18.0 cm	4.0 mm width
30-469-50	18.0 cm	5.0 mm width
30-477-10	23.0 cm	1.0 mm width
30-477-20	23.0 cm	2.0 mm width
30-477-30	23.0 cm	3.0 mm width
30-477-40	23.0 cm	4.0 mm width
30-477-50	23.0 cm	5.0 mm width
30-477-60	23.0 cm	6.0 mm width
30-483-10	20.0 cm	1.0 mm width
30-483-20	20.0 cm	2.0 mm width
30-483-30	20.0 cm	3.0 mm width
30-483-40	20.0 cm	4.0 mm width
30-483-50	20.0 cm	5.0 mm width
30-483-60	20.0 cm	6.0 mm width




Punches

Antrum punch

Working Length: 10.0 cm • Working End: 2.5 mm width
irrigation channel • backwards cutting


Article No.	Working End
30-300-01	straight
30-300-02	left
30-300-03	right



Antrum punch

Working Length: 10.0 cm • Working End: 1.5 mm width
irrigation channel • backwards cutting

Article No.	Working End
30-302-01	straight
30-302-02	left
30-302-03	right



Antrum punch

Working Length: 10.0 cm • irrigation channel
backwards cutting • rotatable


Article No.	Working End
30-304-15	1.5 mm width
30-304-25	2.5 mm width



Antrum punch

Working Length: 10.0 cm • irrigation channel • sideways cutting

Article No.	Working End
30-308-01	left
30-308-02	right



Tobey micro antrum punch

Working Length: 10.0 cm • Shaft: curved • irrigation channel backwards cutting • rotatable

Article No.	Description
30-305-00	downwards
30-307-00	upwards



Hajek-Kofler sphenoid punch

Working Length: 17.0 cm • Working End: 3.0 mm width totally detachable • rotatable

Article No.: 30-456-32




Punches

Hajek-Kofler sphenoid punch

Working Length: 14.0 cm • Working End: 4.0 mm width
totally detachable


Article No.	Working End
30-490-40	90° curved upwards
30-491-40	90° curved downwards



McKenty sphenoid punch

Working Length: 14.5 cm • Working End: 2.0 mm width
throughcutting • rotatable


Article No.: 30-454-20



Sphenoid sinus circular punch

Working Length: 16.0 cm • irrigation channel • 0°


Article No.	Working End
30-422-35	3.5 mm width
30-422-45	4.5 mm width



Sphenoid sinus circular punch

Working Length: 16.0 cm • Working End: straight, egg-form irrigation channel • 60° cutting

Article No.	Working End
30-429-35	3.5 mm width
30-429-45	4.5 mm width



Antrum punch

with collector

Article No.	WL	Working End
30-450-20	18.0 cm	2.0 mm width
30-450-30	15.0 cm	3.0 mm width



Frontal sinus circular punch

Working Length: 16.0 cm • Working End: 60°, egg-form irrigation channel • 90° cut

Article No.	Working End
30-428-35	3.5 mm width
30-428-45	4.5 mm width




Punches

Frontal sinus circular punch

Working Length: 16.0 cm • Working End: 60° • irrigation channel

Article No.	Working End
30-421-35	3.5 mm width
30-421-45	4.5 mm width



Spelemann micro punch

Working Length: 10.0 cm • Shaft: straight

Working End: 2.0 mm width • irrigation channel
backwards cutting • rotatable

Article No.: 30-306-00



Frontal sinus punch

Working Length: 17.0 cm • Working End: 70° curved upwards

Article No.	Working End
30-400-00	3.0 mm width
30-400-01	2.0 mm width



Frontal sinus punch

Working Length: 13.0 cm • Working End: curved upwards, blunt end • irrigation channel

Article No.	Working End
30-425-35	3.5 mm width
30-425-45	4.5 mm width

**Stylet for removal of bone residue**

Working Length: 18.0 cm • Ø: standard (ca, 1.6 mm)
from antrum punch

Article No.: 30-450-90



Micro Instruments

Enucleator



Total Length: 27.0 cm • Shaft: bayonet
Working End: right curved • knurled handel
Article No.: 90-204-02



Hook

Total Length: 27.0 cm • Shaft: bayonet
Working End: curved upwards

Article No.		OD	Working End
90-210-01			blunt
90-210-02	with ball end	1.0 mm	



Hatmann nasal forceps


Article No.	TL	Description
20-200-16	18.0 cm	delicate
20-200-20	20.0 cm	



Krause nasal polyp snare

Total Length: 26.0 cm • Working End: sharp

Article No.: 20-190-25



Nasal snare wires

Working Length: 10.0 cm • OD: 0.35 mm • 100 pieces / box

Article No.: 20-195-35



Cotton applicator

7.5 cm • 5 pieces / box

Article No.: 30-120-07




Cotton carrier

Total Length: 16.0 cm • Working End: 1.2 mm width

5 pieces / box

Article No.	Working End
20-800-16	triangular, serrated
20-801-16	round, serrated
20-802-16	triangular, smooth



Miscellaneous

Injection needle

Total Length: 13.0 cm • Working Length: 10.5 cm • OD: 0.8 mm
Shaft: angled • Luer-Lock

Article No.: 30-122-08



Wigand frontal sinus rasp

Total Length: 21.0 cm • strong


Article No.: 30-209-01



Mirror handle

Total Length: 15.5 cm • Shaft: bayonet • round knurled


Article No.: 90-220-00



Mirror

for use with handle 90-220-00 • Total Length: 12.5 cm
Working End: round • extra delicate


Article No.	OD
90-222-04	4.0 mm
90-222-05	5.0 mm
90-222-06	6.0 mm
90-222-08	8.0 mm
90-222-10	10.0 mm
90-222-12	12.0 mm



Knot tier

Total Length: 20.5 cm • Working End: fork shaped, blunt

Article No.: 40-370-21



Orbita protector

Total Length: 23.0 cm

Article No.: 30-812-24



Key

for opening the slider shaft

Article No.: 20-214-00




Alphabetical Table of Contents

Antrum curette	24	Ritter frontal sinus bougie	12
Antrum curette	24	Rotation clamp	9
Antrum polyp grasping forceps	47	Selfholding Spekulum	8
Antrum punch	54–55	Separator	17
Antrum punch	57	Sickle knife	18
Ball-end probe	13	Sinus frontal curette	27
Binner polyp forceps	46	Sinus speculum	8
Biopsy forceps	49	Sinus suction tube	14
Bipolar forceps	51	Spelemann micro punch	58
Blakesley nasal forceps	33–41, 43	Sphenoid sinus circular punch	56–57
Cottle alar retractor	10	Spoon curette	25–26
Cottle chisel	22	Struempel-Voss nasal forceps	42
Cottle columella forceps	51	Stylet for removal of bone residue	59
Cottle mallet	23	Suction- / Irrigation tube	14
Cottle septum knife	17	Suction- and irrigation sheath	14
Cotton applicator	61	Takahashi nasal forceps	35, 42–43
Cotton carrier	61	Through-cutting nasal forceps	40
Cup forceps	48	Through-cutting suction forceps	41
Dissector	28	Tobey micro antrum punch	55
Draf ethmoid forceps	49–50	Tönnis raspatory	20
Draf ethmoid scissors	31	Uffenorde antrum curette	27
Draf flap knife	17	Walter scissors	30
Draf micro raspatory	20	Weber elevator double ended	19
Draf round knife	17	Weber nasal scissors	31
Draf sickle knife	18	West chisel	22
Eicken antrum cannula	14	West gouge	22
Enucleator	12, 60	Wigand frontal sinus rasp	62
Ferguson suction tube	15	Wigand polyp grasping forceps	42
Fisch retractor	11	Wigand transecting forceps	40
Flexible alligator forceps	48		
Flexible biopsy forceps	48		
Flexible holding arm for use	9		
Frazier suction tube	15–16		
Freer elevator double ended	19		
Frontal sinus biopsy forceps	49		
Frontal sinus circular punch	57–58		
Frontal sinus forceps	45–46, 50		
Frontal sinus forceps	46		
Frontal sinus punch	58–59		
Frontal sinus seeker	12		
Frontal sinus suction raspatory	20		
Gourney suction elevator	19		
Grünwald tubular forceps	44–45		
Guillen suction raspatory	21		
Hajek metal mallet	23		
Hajek-Kofler sphenoid punch	55–56		
Halle antrum curette	24		
Halle antrum curette	24		
Hartmann nasal speculum	7		
Hatmann nasal forceps	61		
Herberhold cranial through-cutting forceps	39		
Heymann nasal scissors	29		
Hook	60		
Howarth elevator double ended	19		
Injection needle	62		
Kerrison sphenoid punch	52–53		
Key	63		
Killian-Struycken septum speculum	7		
Knot tier	63		
Kocher-Langenbeck retractor	10		
Krause nasal polyp snare	61		
Langenbeck-Green retractor	10		
Masing elevator/raspatory	19		
Masing pointed knife	17		
Maxillary sinus seeker	13		
McKenty sphenoid punch	56		
Micro hook	11		
Middeldorpf retractor	10		
Mirror	62		
Mirror handle	62		
Mucosa dissector	28		
Nasal concha scissors	30, 32		
Nasal probe	12		
Nasal snare wires	61		
Orbita protector	63		
Probe hook	13		
Raspatory	20		
Ring curette	27		

10-468-17	11	20-690-03	32	30-200-03	18	30-370-25	44	30-478-10	52	30-828-16	49
20-100-01	7	20-691-01	31	30-200-04	18	30-370-30	44	30-478-20	52	30-828-19	49
20-100-02	7	20-693-00	30	30-203-01	24	30-371-15	44	30-478-30	52	30-832-13	31
20-150-40	7	20-693-01	30	30-203-02	24	30-371-20	44	30-478-40	52	30-832-16	31
20-150-55	7	20-693-02	30	30-205-01	24	30-371-25	44	30-478-50	52	30-832-19	31
20-150-65	7	20-700-10	51	30-205-02	24	30-371-30	44	30-478-60	52	30-833-16	31
20-150-75	7	20-800-16	61	30-205-03	24	30-372-15	44	30-479-10	53	30-834-16	31
20-150-90	7	20-801-16	61	30-205-04	27	30-372-20	44	30-479-20	53	30-835-19	31
20-190-25	61	20-802-16	61	30-205-06	27	30-372-25	44	30-479-30	53	30-943-01	34
20-195-35	61	20-915-00	33	30-205-70	27	30-373-15	45	30-479-40	53	30-943-02	34
20-200-16	61	20-915-01	33	30-207-01	12	30-373-20	45	30-479-50	53	30-943-03	34
20-200-20	61	20-915-02	33	30-207-02	13	30-373-25	45	30-479-60	53	30-943-04	34
20-214-00	63	20-915-03	33	30-207-03	12	30-374-20	45	30-483-10	53	30-943-11	35
20-215-00	36	20-915-04	33	30-209-01	62	30-375-20	45	30-483-20	53	30-943-12	35
20-215-01	36	20-915-10	33	30-300-01	54	30-376-20	45	30-483-30	53	30-943-13	35
20-215-02	36	20-915-12	33	30-300-02	54	30-400-00	58	30-483-40	53	40-370-21	63
20-215-03	36	20-925-00	34	30-300-03	54	30-400-01	58	30-483-50	53	80-126-23	51
20-215-04	36	20-925-01	34	30-302-01	54	30-401-01	46	30-483-60	53	80-306-17	10
20-215-10	36	20-925-02	34	30-302-02	54	30-401-02	46	30-484-10	52	80-306-22	10
20-215-12	36	20-925-03	34	30-302-03	54	30-402-01	46	30-484-20	52	80-306-30	10
20-215-14	36	20-925-04	34	30-304-15	54	30-402-02	46	30-484-30	52	80-308-07	10
20-225-00	37	20-925-11	34	30-304-25	54	30-404-01	46	30-484-40	52	80-308-10	10
20-225-01	37	20-935-01	34	30-305-00	55	30-404-02	45	30-484-50	52	80-308-15	10
20-225-02	37	20-935-02	34	30-306-00	58	30-405-01	46	30-484-60	52	80-308-20	10
20-225-03	37	20-945-01	35	30-307-00	55	30-405-02	50	30-485-10	53	80-334-16	10
20-225-04	37	20-945-02	35	30-308-01	55	30-405-11	50	30-485-20	53	90-201-01	25
20-225-11	37	20-946-01	35	30-308-02	55	30-405-12	50	30-485-30	53	90-201-02	25
20-225-12	37	20-946-02	35	30-310-01	47	30-407-01	50	30-485-40	53	90-201-03	25
20-225-14	37	30-080-05	15	30-310-02	47	30-407-02	50	30-485-50	53	90-201-04	25
20-235-01	37	30-080-06	15	30-311-01	47	30-408-01	50	30-485-60	53	90-201-05	25
20-235-02	37	30-080-07	15	30-311-02	47	30-408-02	50	30-490-40	56	90-201-06	25
20-235-03	37	30-080-08	15	30-320-01	40	30-410-01	49	30-491-40	56	90-201-07	25
20-236-00	42	30-080-09	15	30-320-02	40	30-410-02	49	30-700-00	9	90-201-08	25
20-238-00	42	30-080-10	15	30-321-01	40	30-411-01	49	30-700-01	9	90-201-09	25
20-240-01	41	30-080-12	15	30-321-02	40	30-411-02	49	30-702-25	8	90-201-14	26
20-240-02	41	30-080-15	15	30-322-01	40	30-415-01	46	30-703-25	8	90-201-15	26
20-240-11	41	30-081-06	15	30-322-02	40	30-415-02	46	30-705-55	8	90-201-16	26
20-240-12	41	30-081-08	15	30-328-01	48	30-421-35	58	30-705-65	8	90-201-17	26
20-245-02	42	30-081-09	15	30-330-01	40	30-421-45	58	30-705-75	8	90-201-18	26
20-245-11	42	30-081-10	15	30-330-02	40	30-422-35	56	30-706-55	8	90-201-19	26
20-245-12	42	30-081-12	15	30-332-01	39	30-422-45	56	30-706-65	8	90-201-28	26
20-246-02	42	30-081-15	15	30-332-02	39	30-425-35	59	30-706-75	8	90-201-29	26
20-300-22	19	30-082-08	15	30-333-01	39	30-425-45	59	30-711-00	8	90-202-01	18
20-305-18	19	30-082-10	15	30-333-02	39	30-428-35	57	30-740-10	20	90-202-02	18
20-306-20	19	30-082-12	15	30-340-01	38	30-428-45	57	30-740-30	20	90-204-01	12
20-307-18	19	30-082-15	15	30-340-02	38	30-429-35	57	30-740-31	20	90-204-02	60
20-313-00	12	30-084-01	16	30-340-03	38	30-429-45	57	30-740-40	20	90-206-01	20
20-315-21	19	30-084-02	16	30-342-01	38	30-450-20	57	30-742-19	18	90-206-02	20
20-320-19	19	30-100-25	14	30-342-02	38	30-450-30	57	30-743-19	17	90-206-03	20
20-322-01	21	30-100-30	14	30-342-03	38	30-450-90	59	30-744-19	17	90-206-04	20
20-322-02	21	30-100-40	14	30-343-01	38	30-454-20	56	30-747-06	15	90-208-01	13
20-322-03	21	30-101-25	14	30-343-02	38	30-456-32	55	30-747-07	15	90-208-02	13
20-323-01	21	30-101-30	14	30-343-03	38	30-468-10	52	30-747-08	16	90-210-01	60
20-323-02	21	30-101-40	14	30-343-04	38	30-468-20	52	30-747-09	16	90-210-02	60
20-323-03	21	30-102-25	14	30-343-11	39	30-468-30	52	30-747-10	16	90-210-03	28
20-400-18	23	30-102-30	14	30-343-12	39	30-468-40	52	30-747-12	16	90-210-04	27
20-402-21	23	30-102-40	14	30-343-13	39	30-468-50	52	30-747-15	16	90-210-05	27
20-405-03	22	30-103-15	16	30-343-14	39	30-469-10	53	30-749-15	14	90-210-06	27
20-405-04	22	30-105-06	16	30-344-00	41	30-469-11	52	30-759-21	20	90-210-07	27
20-405-06	22	30-105-08	16	30-344-45	41	30-469-20	53	30-805-00	24	90-210-08	27
20-405-07	22	30-105-10	16	30-346-00	43	30-469-21	52	30-807-01	24	90-210-09	17
20-405-09	22	30-106-25	14	30-346-30	43	30-469-30	53	30-807-02	24	90-210-10	17
20-405-10	22	30-106-30	14	30-350-00	43	30-469-31	52	30-807-03	24	90-211-01	26
20-405-12	22	30-106-40	14	30-350-01	43	30-469-40	53	30-809-01	27	90-212-01	26
20-500-13	17	30-108-00	20	30-350-02	43	30-469-41	52	30-809-02	27	90-212-02	26
20-501-14	17	30-120-07	61	30-350-03	43	30-469-50	53	30-812-24	63	90-212-03	26
20-502-14	17	30-121-25	12	30-352-00	43	30-469-51	52	30-816-40	22	90-213-00	13
20-560-14	10	30-121-30	12	30-352-01	43	30-471-10	53	30-817-40	22	90-214-00	11
20-560-15	10	30-121-40	12	30-352-02	43	30-471-20	53	30-820-10	50	90-215-00	28
20-561-15	10	30-121-50	12	30-354-01	43	30-471-30	53	30-820-13	50	90-220-00	62
20-630-10	30	30-122-08	62	30-354-02	43	30-473-10	52	30-820-16	50	90-222-04	62
20-650-16	29	30-125-00	49	30-358-00	40	30-473-20	52	30-820-19	50	90-222-05	62
20-651-18	29	30-164-00	14	30-358-01	40	30-473-30	52	30-822-13	50	90-222-06	62
20-680-16	29	30-167-00	14	30-358-02	40	30-477-10	53	30-822-16	50	90-222-08	62
20-681-16	29	30-167-30	14	30-360-00	40	30-477-20	53	30-822-19	50	90-222-10	62
20-682-16	29	30-172-40	48	30-360-01	40	30-477-30	53	30-824-13	50	90-222-12	62
20-683-18	29	30-173-40	48	30-360-02	40	30-477-40	53	30-824-16	50		
20-690-01	32	30-200-01	18	30-370-15	44	30-477-50	53	30-824-19	50		
20-690-02	32	30-200-02	18	30-370-20	44	30-477-60	53	30-826-16	50		











Charrière (CH) French (F)	mm	Inch
1	0.33	0.013
2	0.67	0.025
3	1.00	0.039
4	1.33	0.051
5	1.67	0.065
6	2.00	0.079
7	2.33	0.092
8	2.67	0.105
9	3.00	0.118
10	3.33	0.131
11	3.67	0.144
12	4.00	0.156
13	4.33	0.170
14	4.67	0.184
15	5.00	0.197
16	5.33	0.210
17	5.67	0.223
18	6.00	0.237
19	6.33	0.249
20	6.67	0.263
21	7.00	0.276
22	7.33	0.288
23	7.67	0.302
24	8.00	0.315
25	8.33	0.328
26	8.67	0.341
27	9.00	0.354
28	9.33	0.367
29	9.67	0.380
31	10.23	0.407
32	10.67	0.420
34	11.33	0.446
36	12.00	0.472
38	12.67	0.494
40	13.33	0.520
45	15.00	0.591



For products manufactured by SPIGGLE & THEIS Medizintechnik GmbH, the electronic instructions for use are available at the following website:

stm-ifu.com

You will need a PDF reader to view the documents. If you do not have one installed on your device, you can download a free reader, for example from Adobe. If a paper version is required, SPIGGLE & THEIS Medizintechnik GmbH will provide a printed copy free of charge within seven days of receiving the request. Please contact us at the following email address: ifu@spiggle-theis.com.

